

# Astrella

### Ultrafast Ti:Sapphire Amplifier

Astrella and the new Astrella HE are next-generation, ultrafast, kHz amplifiers that are the first to combine industry-leading performance and industrialized durability. Manufactured to Coherent's rigorous standards using advanced stress-testing techniques, the one-box Astrella system enables a wide range of demanding scientific applications and operating conditions, offering higher productivity and lower data acquisition costs. Delivering high (up to >9 mJ/pulse) energy, either <35 fs or <100 fs pulse widths, and excellent beam quality (M² <1.25), Astrella is ideal for ultrafast spectroscopy, THz studies, femtosecond micromachining, etc. With unmatched performance, reliability and affordability, Astrella stands at the forefront of the industrial revolution in ultrafast science.



#### **FEATURES & BENEFITS**

- · One-box, industrialized platform
- HASS\* verified for quality and reliability
- >5 mJ, >7 mJ or >9 mJ <35 fs or <100 fs pulses</li>
- High performance STAR regen amplifier (water-only cooling)
- · Hands-free Vitara oscillator
- Revolution pump laser for performance overhead
- Sealed stretcher/compressor section with advanced dispersion management for clean, short pulses
- Thermally-stabilized sub-systems for long term stability

#### **APPLICATIONS**

- Time-resolved Spectroscopy
- Multidimensional Spectroscopy
- THz Spectroscopy
- fs Micromachining
- Surface SFG/SHG
- · Stimulated Raman Scattering
- \* HASS Highly Accelerated Stress Screening




SPECIFICATIONS <sup>1</sup>	Astrella USP	Astrella F	Astrella HE USP	Astrella HE F
Center Wavelength <sup>2</sup> (nm) (nominal)	795 to 805	780 to 820	795 to 805	780 to 820
Repetition Rate <sup>3</sup> (kHz)	1,5			
Pulse Duration <sup>3,4</sup> (fs) (FWHM)	<35	<100	<35	<100
Contrast Ratio⁵				
Pre-Pulse	>1000:1			
Post-Pulse	>100:1			
Power Stability <sup>6,7</sup> (rms)	<0.5			
Beam Pointing Stability <sup>6,7</sup> (mrad) (rms)	<10			
Beam Diameter (mm) (1/e <sup>2</sup> ) (nominal)				
1 kHz	11		13	
5 kHz			1	
Spatial Mode	$TEM_{00}, M^2 < 1.25$			
Polarization	linear, horizontal			
Energy per Pulse (mJ)				
1 kHz	>5.0, >7.0		>9.0	
5 kHz	>1	.4	>2.0	
Pump Laser	Revolution-50,	Revolution-65	Revolution-80	
Seed Laser	Vitara-S, Vitara-T, or Vitara-T-HP			
Each System HASS Verified	Yes			


Specifications apply at 800 nm.

#### TYPICAL PERFORMANCE DATA

#### **Astrella HE 24-Hour Stability**



#### **Astrella HE USP Pulse Width**





Factory set, must be specified when ordered and will be optimized prior to shipment.
 Contact factory for other repetition rates and pulse width options.

<sup>4</sup> A Gaussian pulse shape de-convolution factor (0.7) is used to determine the pulse width from an autocorrelation signal measured by a Coherent SSA (Single-Shot Autocorrelator).


<sup>5</sup> Contrast ratio is defined as the ratio between the peak intensity of the output pulse to the peak intensity of any other pulse that occurs greater than 1 ns before or after the output pulse.
6 Under stable environmental conditions after system warm-up.

<sup>7</sup> Over 24 hrs.


### **MECHANICAL SPECIFICATIONS**

#### Astrella


#### **Right Side View**





#### **Top View**



## Rear View









Coherent, Inc., 5100 Patrick Henry Drive Santa Clara, CA 95054 p. (800) 527-3786 | (408) 764-4983 f. (408) 764-4646

tech.sales@coherent.com www.coherent.com