

Chameleon Compact OPO

Automated Wavelength Extension for Chameleon Ti:Sapphire Lasers




Features

- Automated hands-free wavelength tuning
- Independent wavelength tuning of both laser and OPO
- Two output ports for Ti:Sapphire and OPO
- Fan poled OPO crystal for flexible pump wavelength tuning
- Full tuning range without change of optics or re-alignment


- Gap-free overlap with edge of Ultra II tuning range
- USB and RS-232 control
- Optimized for pumping with the Chameleon Ti:Sapphire family
- Effortless tuning over water absorption bands
- Idler option available
- Remote service facility

Mechanical Specifications

Chameleon Compact OPO - Top View



Chameleon - Side View




Chameleon™ Compact OPO

Automated Wavelength Extension for Chameleon Ti:Sapphire Lasers


System Specifications		Pumped by Chameleon Ultra II	Pumped by Chameleon Vision II	
	Tuning Range¹ (nm)	1000 to 1600	1000 to 1350	
	Pump Wavelength Range (nm)	740 to 880	760 to 880	
	Output Power, Signal ² (mW)	>700	>700³	
	Pump Output Power Available ⁴ (%) when pumping OPO in bypass mode	when pumping OPO ~15 n bypass mode ~95 width (fs)(typical) 200 /pical) 1.1 Diameter (mm) 2		
	Pulsewidth (fs)(typical)			
	M² (typical)			
	Beam Diameter (mm)			
	Beam Divergence (mrad)(typical)			
	Polarization	Horizontal		
	Repetition Rate (MHz)	8o, locked to pump laser		
	Dimensions (L x W x H)	520 x 369 x 158 mm (20.5 x 14.5 x 6.2 in.)		
Idler Option	Additional Tuning Range ^{1,5} (nm)	1750 to 4000		
	Idler Output Power ⁶ (mW)	>100		
Max IR Option	Additional Signal Output Power ⁷ (%)	1	5	

- ¹ Tuning range depends on Pump Wavelength.
- ² At maximum of pump and OPO signal tuning curve.
- 3 Vision Dispersion settings optimized for OPO.
- 4 Typical. Please refer to Chameleon datasheet for respective power specifications. Pump output not available with Max IR option.
- ⁵ In addition to OPO signal output range. All other specifications are unaffected.
- ⁶ At maximum of pump and OPO idler tuning range .
- ⁷ No access to pump beam.

Chameleon Ultra II and Compact OPO Tuning Range



Chameleon Compact OPO Typical Output Power



Coherent follows a policy of continuous product improvement. Specifications are subject to change without notice

Coherent's scientific and industrial lasers are certified to comply with the Federal Regulations (21 CFR Subchapter J) as administered by the Center for Devices and Radiological Health on all systems ordered for shipment after August 2, 1976.

Coherent offers a limited warranty for all Chameleon systems. For full details of this warranty coverage, please refer to the Service section at www.Coherent.com or contact your local Sales or Service Representative.



Coherent, Inc. 5100 Patrick Henry Drive Santa Clara, CA 95054 phone (800) 527-3786

phone (800) 527-3786 (408) 764-4983

fax (408) 764-4646 e-mail tech.sales@Coherent.com Benelux +31 (30) 280 6060 China +86 (10) 8215 3600 France +33 (0)1 8038 1000 Germany +49 (6071) 968 333 Italy +39 (02) 31 03 951

Japan +81 (3) 5635 8700 Korea +82 (2) 460 7900 UK +44 (1353) 658 833

