
Q-SMART 1200 & 1500

Compact High-Energy pulsed Nd:YAG lasers with excellent beam quality and versatility





MAIN FEATURES


- Up to 1.5 J @ 1064 nm
- Robust and field proven technology
- Design to last thanks to ceramic reflectors and 100 million shots flashlamp lifetime warranty
- Plug & play harmonic modules with automatic phase-matching
- Cables and cooling lines fully disconnectable
- · Easy to use and maintain
- · No need for external water
- Universal voltage
- Intuitive GUI interface
- · SLM option (Single Longitudinal Mode)

MAIN APPLICATIONS

- LiDAR
- INSTRUMENTATION
- · PLI
- · DYE, OPO AND TI:SA PUMPING
- SPECTROSCOPY
- LIF AND COMBUSTION

•••


Typical beam profiles



Near Fied 1.5 J @ 1064 nm, 10 Hz



Far Field 1.5 J 1064 nm, 10 Hz



Temporal Profile @ 1064 nm (1 GHz scope)

www.quantel-laser.com

Please contact Lumibird to find the best match fo your needs and compatibility between options.









Q-SMART 1200 & 1500

Compact High-Energy pulsed Nd:YAG lasers with excellent beam quality and versatility



SPECIFICATIONS

		Q-smart 1200	Q-smart 1500
Repetition rate (1) (Hz)		10	10
Energy per pulse (mJ)	1064 nm	1200	1500
	532 nm ⁽²⁾	600	820
	355 nm ⁽²⁾	280	490
	266 nm	130	150
Pulse duration(3) (ns)	1064 nm	5-10	
Beam diameter (mm)	1064 nm	≤ 10	
Divergence ⁽⁴⁾ (mrad)	1064 nm	≤ 0.5	
Polarization(5) ratio (%)	1064 nm	≥ 80	
Spatial profile	Near field ⁽⁶⁾ (fit to Gaussian)	≥ 0.7	
	Far field ⁽⁷⁾ (fit to Gaussian)	≥ 0.9	

Service requirements		
Power	2 x 100-240 VAC 50-60 Hz Single phase	
Cooling group	Air to water	


Laser head

- Other repetition rates on request
- (2) HE: High Energy option on request


- (3) Measured at FWHM with fast photodiode and 1 GHz scope
 (4) Full angle at 1/e² of peak
 (5) Horizontal @ 1064 nm Vertical @ 532 nm Horizontal @ 355 nm and @ 266 nm
- (6) Measured at 1 m from laser output
- (7) Measured at focal plane of a 2 m focus lens, least square fit to Gaussian (perfect fit = 1)

Power drift ⁽¹⁾ (%)	1064 nm	± 3
	532 nm	± 5
	355 nm	± 5
	266 nm	± 10
Energy stability ⁽²⁾ (%)	1064 nm	± 2 (0.6)
	532 nm	± 4 (1.3)
	355 nm	± 6 (2)
	266 nm	± 8 (2.6)
Pointing stability(3) (µrad)	1064 nm	< 40
Jitter ⁽⁴⁾ (ns)	Standard	± 0.5
	SLM option	± 1
M², focusability (times diffraction limit)	1064 nm	≤ 2
Linewidth (cm ⁻¹)	Standard ⁽⁵⁾	≤ 0.7
	SLM ⁽⁶⁾ option	≤ 0.005
Temperature range	18-28°C	

- (1) Over 8 hours for $\Delta T^{\circ} \leq \pm 3^{\circ}C$
- (2) Peak-to-peak (RMS), 99% of shots
- (3) Measured with Spiricon LBA-100, rms, on 200 pulses at the focal plane
- of a 1 m focus lens
- (4) With respect to Q-Switch trigger, at half-width of 500 accumulated shots for 99% of shots
- (5) Measured at FWHM with a grating spectrometer with 0.045 cm-1 resolution
- (6) Measured with a slow scan Fabry-Perot Etalon, ≤ 20% energy reduction @ 1064 nm



Electronics



www.quantel-laser.com







