

Harmonic Generator System (HGS)


Second-, Third-, and Fourth- Harmonic Accessory for Ti:Sapphire Amplifiers





Features

- Simple collinear mixing configurations
- Compatible with Coherent 35 fs to 130 fs, and 0.5 ps to 1.5 ps Ti:Sapphire amplifiers (Legend Elite™, Libra™, Hidra)
- High-damage-threshold, high-conversion efficiency BBO crystals
- Includes wavelength separation mirrors, waveplates, focusing optics, and more
- Stand-alone enclosure two enclosures are required for fourth harmonic

Mechanical Specifications







Harmonic Generator System (HGS)

Second-, Third-, and Fourth- Harmonic Accessory for Ti:Sapphire Amplifiers

Specifications		HGS-D	HGS-T	HGS-F	
	Conversion Efficiency (%, from 800 pump, fsec)	35	15	1	
	Power Stability (%, rms)	1.5	2.5	3.5	

Coherent follows a policy of continuous product improvement. Specifications are subject to change without notice.

Coherent's scientific and industrial lasers are certified to comply with the Federal Regulations (21 CFR Subchapter J) as administered by the Center for Devices and Radiological Health on all systems ordered for shipment after August 2, 1976.

Coherent offers a limited warranty for all Harmonic Generator Systems. For full details of this warranty coverage, please refer to the Service section at www.Coherent.com or contact your local Sales or Service Representative.



Coherent, Inc. 5100 Patrick Henry Drive Santa Clara, CA 95054 phone (800) 527-3786 (408) 764-4983 fax (408) 764-4646

fax (408) 764-4646 e-mail tech.sales@Coherent.com Benelux +31 (30) 280 6060
China +86 (10) 8215 3600
France +33 (0)1 8038 1000
Germany +49 (6071) 968 333
Italy +39 (02) 31 03 951
Japan +81 (3) 5635 8700
Korea +82 (2) 460 7900
UK +44 (1353) 658 833