

PRIORITY	
_____ Mandatory	_____ As Required
_____ Next Visit	_____ At Installation
<u>  X  </u> Information	

# FIELD SERVICE BULLETIN

**DATE ISSUED:** April 25, 2014

**NUMBER:** 733


## Coherent GUI v2.41 Error

**DATE ISSUED** 25<sup>th</sup> April 2014

**PRODUCT** Coherent GUI v2.41 32bit & 64bit

**PURPOSE** Inform Field Service Engineers of an error with the Coherent GUI v2.41 when used with Chameleon lasers systems. This occurs after the software is downloaded from the Coherent website.

**DESCRIPTION** 'The download of the JRE failed' appears while opening the Coherent GUI v2.41 after downloading the software from the Coherent website. The GUI can still be opened, however the spectrometer will not function and other random issues may occur.



**SOLUTION** Use the 'GLA transition GUI' which can be found within Spectrum <https://spectrum.na.coherentinc.com/service/home/dpssScotlandSVC/SitePages/Home.aspx>.

**For more information contact:** PSE OM