

Chameleon MPX


Multiphoton Microscopy Wavelength Extension

Chameleon MPX extends the wavelength range of Chameleon Vision and Ultra Ti:Sapphire lasers, and is specifically designed and optimized for non-linear imaging techniques.

Employing the latest generation fan-poled OPO technology, the fully automated Chameleon MPX delivers high peak power to the sample plane with short pulses and dispersion compensation optimized for typical commercial microscope systems.

Featuring a wide pump tuning range, the Chameleon MPX offers independently tunable dual beam excitation of popular fluorescent probes (e.g. eGFP, mCherry), enabling powerful and truly flexible multimodal imaging.

**Chameleon MPX
Dual Tuning Range**



Superior Reliability & Performance

Chameleon MPX Features:

- **Performance optimized for non-linear imaging applications**
- **130 fs short pulses for high peak power**
- **GDD dispersion compensated output**
- **Automated, gap-free tuning from 680 nm to 1340 nm with Chameleon pump laser**
- **Independent wavelength tuning of both laser and OPO for simultaneous 2-color excitation and multimodal imaging**
- **Synchronized output pulse trains for CARS/SRS and wavelength mixing**

Chameleon MPX Applications:

- **Multiphoton Excitation Microscopy**
- **Second Harmonic Generation Imaging**
- **Third Harmonic Generation Imaging**
- **CARS/SRS Microscopy**
- **Ultrafast Spectroscopy**
- **Non-linear Optics**

Chameleon MPX

Multiphoton Microscopy Wavelength Extension

System Specifications

	Pumped by Chameleon Ultra II	Pumped by Chameleon Vision II
Tuning Range ¹ (nm)	1010 to 1340	1010 to 1340
Pump Wavelength Range (nm)	760 to 910	760 to 910
Output Power ² (mW)(signal)	>750	>700 ³
Pump Output Power Available ⁴ (%)		
when pumping OPO		~15
in bypass mode		95
Pulse Width ⁵ (fs)(typical)		130
Output GDD Precompensation ⁶ (fs ²)		-6000
M ² (typical)		<1.1
Beam Diameter (mm)		2
Beam Divergence (mrad)(typical)		0.7
Polarization		Horizontal
Repetition Rate (MHz)		80 MHz, locked to pump laser
Dimensions (L x W x H)		520 x 369 x 158 mm (20.5 x 14.5 x 6.2 in.)

¹ Tuning range depends on Pump Wavelength.

² At maximum of pump and OPO signal tuning curve.


³ Vision Dispersion settings optimized.

⁴ Typical. Please refer to Chameleon datasheet for respective power specifications.

⁵ Typical value at sample plane after microscope dispersion.

⁶ Typical value at 1100 nm.

Chameleon MPX Typical Output Power




Chameleon MPX


Multiphoton Microscopy Wavelength Extension

Mechanical Specifications

Chameleon MPX - Top View



Chameleon - Side View



Chameleon and Chameleon MPX Table Layout



COHERENT®

www.Coherent.com

Coherent, Inc.,

5100 Patrick Henry Drive
Santa Clara, CA 95054

phone (800) 527-3786
(408) 764-4983

fax (408) 764-4646

e-mail tech.sales@Coherent.com

Benelux +31 (30) 280 6060

China +86 (10) 8215 3600

France +33 (0)1 8038 1000

Germany/Austria/

Switzerland +49 (6071) 968 333

Italy +39 (02) 31 03 951

Japan +81 (3) 5635 8700

Korea +82 (2) 460 7900


Taiwan +886 (3) 505 2900

UK/Ireland +44 (1353) 658 833

Coherent follows a policy of continuous product improvement. Specifications are subject to change without notice.

Coherent's scientific and industrial lasers are certified to comply with the Federal Regulations (21 CFR Subchapter J) as administered by the Center for Devices and Radiological Health on all systems ordered for shipment after August 2, 1976.

Coherent offers a limited warranty for all Chameleon systems. For full details of this warranty coverage, please refer to the Service section at www.Coherent.com or contact your local Sales or Service Representative.



CE ISO 9001 Registered